



To Whom It May Concern:

Thank you for using HeartMath equipment. Below are the cleaning instructions.

What are the infection control and cleaning instructions for the HeartMath sensors?

Alcohol wipes and the germicidal wipes are suitable for sanitizing the Inner Balance ear sensors, emWave Pro ear and finger sensors, and emWave2 ear sensors.

Be sure that no liquid gets inside the sensor.

Please let us know if you have any questions or would like further information.

Sincerely,

Tom Beckman
HeartMath LLC
831-338-8745
tom@heartmath.com